



Taha Kimia TajhizCo.



Core Lab



Manual Gas 100 Datasheet

Reservoir Fluid Analysis-PVT \ Gasometry



Manual Gas 100



Specifications	
Volume	100 ml
Pressure	Atmospheric
Temperature	Ambient
Volume Accuracy	± 0.1 ml
Pressure Accuracy	± 0.1 mbar
Temperature Accuracy	$\pm 0.1^\circ$

General Features

This equipment is composed of following elements :

- Absolute Linear Coder
- 2 Manual Valves
- Absolute Pressure sensors

All wetted parts are made of Stainless Steel (ref SS316).

