



Taha Kimia TajhizCo.



Core Lab



# Volumetric Manual Pump VMPSS 30-700 Datasheet

Pump\ Special Pumps



## Volumetric Manual Pump Single Screw – VMPSS 30-700



Specifications	
Volume	30 ml
Pressure	700 bar
Temperature	Ambient
1 <sup>st</sup> Gear	0.6 ml/rev
2 <sup>nd</sup> Gear	0.02 ml/rev

### General Features

How to inject manually under high pressure with no effort?

The manual syringes pump VMPSS 30-700 is specially designed to inject or control pressure manually with no effort. The two stage gear box permits the injection at high flow in low pressure and decrease or control at high pressure with the special gear box by manual switching between high and low speed.

The concept permits fillings of pump or system in short time and control, with a good accuracy and no effort, the pressure. The system is used for systems testing, filling as well cleaning.