



Taha Kimia TajhizCo.




Core Lab



# Manual GOR System Datasheet

Reservoir Fluid Analysis-PVT\Gas Oil Ratio

# GOR System



GOR System is dedicated to measure gas/oil ratio by flashing pressurized liquids at equilibrium conditions. This system is composed of a panel supporting a graduated glass tube, a density cell, a liquid trap cell and a sampling bottle. This easy to use GOR system ensures accurate measurement of liberated gas thanks to the inlet/outlet connection for gasmeter and gravity technique allows measurement of flashed dead oil.

Specifications	
Working Pressure	Atm
Working Temperature	Ambient
Fluids	Hydrocarbons G1