



Taha Kimia TajhizCo.



Core Lab



QUAD Viscometer Datasheet

Reservoir Fluid Analysis-PVT\Onfield Instruments



QUAD Viscometer



General Features

The obligation to obtain data on field as fast as possible for its exploitation led us to realize in collaboration with and at the request of petroleum companies a transportable unit for on-site PVT analysis.

With QUAD equipment, proposed in the On-field Instruments section, we propose a viscometer.

Specifications	
Working Pressure	1000 bar
Range	0.1 to 10 000 cp
Working Temperature	Ambient to 200°C
Volume	10 ml
Size	1020 x 420 x 380 mm



QUAD Viscometer

Why choosing ST QUAD Viscometer?

- Capillary Viscometer
- Fast Connections
- Determination of Viscosity
- Data Acquisition and Full Automation with Falcon® Software
- Automatic Control of Pressure / Flow / Temperature

