



Taha Kimia Tajhiz Co.



Core Lab



Automated Acid Attack - AAA Datasheet

Geochemistry



Automated Acid Attack - AAA



General Features

This equipment has been made to realize samples demineralization thanks to acid attacks. All the process is automated and totally secure for the user. Preparation as Kerogen or Palynologic can be prepared with this AAA.

Specifications	
Working Temperature	+20°C to 80°C
Number of Reactor	8
Reactor Volume	200 and 500 ml
Sample Volume	Up to 65 ml



Automated Acid Attack - AAA

Our AAA is composed of following parts:

- 3 Stigma pump: 1 for Acid Injection, 1 for Water injection and 1 for collection
- 8 Reactor in staggered rows for an easy access
- Climatic Air Bath with 2 door glass.
- Effluent Recovery System
- Falcon® Software

Example of Synoptic

