



Taha Kimia Tajhiz Co.



OFITE

# BD7 Trim Saw

## Datasheet

**Core Preparation / Slabbing-Trimming-Grinding**




## BD7 Trim Saw

- Large 14" x 14" cutting surface
- Direct drive motor
- All cast aluminum tank, table and blade guard
- Table tilts up and out of the way for blade removal/change, adding oil, or cleaning
- Easy to clean tank with oil drain hose and sludge removal plug
- Mounted on a 1/2" polypropylene base
- 7" Professional Lapidary Diamond Blade

### Specifications

Horse Power	1/3
Voltage*	115V/60Hz
Amperage	4.4
Motor RPM	3,300
Blade Arbor	5/8"
Blade Capacity	7"
Depth of Cut	1-3/16"
Max. Vise Opening	2-3/4"
L x W x H (in)	17" x 18" x 12"
L x W x H (mm)	432 x 457 x 305 mm
Saw Weight	28 lbs.
<i>With vise &amp; splash hood</i> Part #	168765
<i>Without vise &amp; splash hood</i> Part #	170789

\*220Volt Part # 169780 (with vise & splash hood)



Built in splash guard



Optional:  
PETG high-impact,  
copolymer splash hood



Optional:  
Adjustable rock vise  
clamp system cross  
feed adjustment

For more information  
800.630.7682