



Taha Kimia TajhizCo.



Core Lab



Standard Sample for Permeability Calibration SSPC Datasheet

Core Preparation/ Calibration (Check Plugs)



Standard Sample for Permeability Calibration



General Features

For checking and testing the calibration of gas permeameters and porosimeters, St uses Standard Sample for Permeability Calibration. This set is composed of three plugs with a different permeability and porosity range.

St Standard Samples are based on real rock samples whose permeability and porosity physical specifications are given through testing report from a geology laboratory. These detailed reports are included in the set.

Specifications	
Diameter	1½"
Length	2"
Permeability Range	0.05 mD to 250 mDa
Porosity Range	From 4 to 15%