



Taha Kimia Tajhiz Co.



Core Lab


# Contact Angle Interfacial Tension, IFT-10-C Datasheet

Fluid Characteristics

#



## Contact Angle Interfacial Tension IFT-10-C



A contact angle measurement device can be added to the Pendant Drop system IFT 10-P. The contact angle measurement is made by placing a liquid drop on a piece of quartz, core, or other test material. The angle of contact between the drop and the test piece is measured under reservoir pressure and temperature conditions. This system can also be supplied as a stand alone system.